



Recombinant African swine fever virus Virus attachment protein p12 (Pret-110)-VLPs

Product Code	CSB-MP315474AYG
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Uniprot No.	P0C9Y3
Form	Lyophilized powder
Storage Buffer	Lyophilized from PBS, 6% Trehalose, pH 7.4
Product Type	Transmembrane Proteins
Immunogen Species	African swine fever virus (isolate Tick/South Africa/Pretoriuskop Pr4/1996) (ASFV)
Sequence	MALDGSSGGGSNVETLLIVAIIVVIMAIMLYYFWWMPRQQKKCSKAEECTCNN GSCSLKTS
Research Area	Others
Source	Mammalian cell
Target Names	Pret-110
Protein Names	Recombinant African swine fever virus Virus attachment protein p12(Pret-110)
Expression Region	1-61aa
Notes	Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity.
Tag Info	C-terminal 10xHis-tagged If you have specified tag type, please tell us and we will check if it's possible to develop.
Length/Coupling Ratio	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.