



Recombinant Sceloporus woodi L-lactate dehydrogenase B chain (LDHB)

Product Code	CSB-EP012841SDH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79913
Product Type	Recombinant Protein
Immunogen Species	Sceloporus woodi (Florida scrub lizard)
Purity	>85% (SDS-PAGE)
Sequence	ASLKDKLIT PVAQPATQPT SKVTVVGVGQ VGMACAISVL EKGLCDELAL VDVLEDKLLKG EMMDLQHGS LFLKTNKIVAG KDYAVTANSK VVVVTVAGVRQ QEGESRLDLV QRNVNVFKFI IPQVVKYSPD CIILVVSNPV DILTYVTWKL SGLPKHRVIG SGCNLD SARF RFLMGERLGI HPSSCHGWIL GEHGDSSVAV WSGVNVAGVS LQELNPAMGS DQDSEGWKQV HKQVVD SAYE VIKLLKGYTNW AIGLSVADLL ETIMKNLCRV HPVSTMVKGM YGIENEVFLS LPCVLG SVGL TSVINQKLKD SEVAQLQ TSA TTLWNVQKDL KDL
Source	E.coli
Target Names	LDHB
Protein Names	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
Expression Region	2-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.