



# Recombinant Sceloporus woodi L-lactate dehydrogenase B chain (LDHB)

<b>Product Code</b>	CSB-EP012841SDH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P79913
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sceloporus woodi (Florida scrub lizard)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASLKDKLIT PVAQPATQPT SKVTVVGVGQ VGMACAISVL EKGLCDELAL VDVLEDKCLKG EMMDLQHGS LFLKTNKIVAG KDYAVTANSK VVVVTVAGVRQ QEGESRLDLV QRNVNVFKFI IPQVVKYSPD CIILVVSNPV DILTYVTWKL SGLPKHRVIG SGCNLSARF RFLMGERLGI HPSSCHGWIL GEHGDSSVAV WSGVNVAGVS LQELNPAMGS DQDSEGWKQV HKQVVD SAYE VIKCLKGYTNW AIGLSVADLL ETIMKNLCRV HPVSTMVKGM YGIENEVFLS LPCVLGVSGL TSVINQKLKD SEVAQLQ TSA TTLWNVQKDL KDL
<b>Source</b>	E.coli
<b>Target Names</b>	LDHB
<b>Protein Names</b>	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
<b>Expression Region</b>	2-333
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.