




Recombinant Human Glucagon receptor (GCGR)-VLPs

Product Code	CSB-MP009316HU
Abbreviation	Recombinant Human GCGR protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47871
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	The purity information is not available for VLPs proteins.
Sequence	AQVMDFLFEKWKLYGDQCHHNSLLPPPELVCNRTFDKYSCWPDTPANTTA NISCPWYLPWHHKVQHRFVFKRCGPDGQWVRGPRGQPWRDASQCQMDGE EIEVQKEVAKMYSSFQVMYTVGYSLSLGALLLAILGGLSKLHCTRNAIHANLF ASFVLKASSVLVIDGLLRTRYSQKIGDDLSTWLSGDGAVAGCRVAAVFMQYG IVANYCWLLVEGLYLHNLGLATLPERSFFSLYLGIGWGAPMLFVVPWAVVKCL FENVQCWTSNDNMGFWWILRFPVFLAILINFFIFVRIVQLLVAKLRARQMHTD YKFR LAKSTLTIPLLVHEVVFAFVTDEHAQGTLRS AKLFFDLFLSSFQGLLVA VLYCFLNKEVQSELRRRWHRWRLGKVLWEERNTSNHRASSSPGHGPPSKEL QFGRGGGSQDSSAETPLAGGLPRLAESPF
Research Area	Cancer
Source	Mammalian cell
Target Names	GCGR
Expression Region	26-477aa
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	52.6 kDa
Protein Length	Full Length of Mature Protein
Image	 CSB-MP009316HU is detected by Mouse anti-6*His monoclonal antibody
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize



for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.