



Recombinant Mycoplasma genitalium Uncharacterized protein MG281 (MG281)-VLPs

Product Code	CSB-MP341708MLN
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Uniprot No.	P47523
Form	Lyophilized powder
Storage Buffer	Lyophilized from PBS, 6% Trehalose, pH 7.4
Product Type	Transmembrane Proteins
Immunogen Species	Mycoplasma genitalium (strain ATCC 33530 / G-37 / NCTC 10195)
Sequence	<p>MQFKKHKNSVKFKRKLFWTIGVLGAGALTTFSAVMITNLVNQSGYALVASGRS GNLGFKLFSTQSPSAEVKLKSLSLNDGSYQSEIDLSSGGANFREKFRNFANELS EAITNSPKGLDRPVPKTEISGLIKTGDNFITPSFKAGYYDHVASDGSLLSYYQST EYFNNRVLPILQTTNGTLMANNRGYDDVFRQVPSFSGWSNTKATTVSTSNN LTYDKWTFYFAAKGSPLYDSYPNHFFEDVKTLAIDAKDISALKTTIDSEKPTYLIIR GLSGNGSQLNELQLPESVKKVSLYGDYTGNNVAKQIFANVVELEFYSTSKANS FGFNPLVLGSKTNVIYDLFASKPFTHTIDLTQVTLQNSDNSAIDANKLKQAVGDIY NYRRFERQFQGYFAGGYIDKYLKVNVTNKDSDDDLVIYRSLKELNLHLEEAYR EGDNTYYRVNENYYPGASIIYENERASRDSEFQNEILKRAEQNGVTFDENIKRIT ASGKYSVQFQKLENDTDSSLERMTKAVEGLVTVIGEEKFETVDITGVSSDTNE VKSLAKELKTNALGVKLLK</p>
Research Area	Others
Source	Mammalian cell
Target Names	MG281
Protein Names	Recombinant Mycoplasma genitalium Uncharacterized protein MG281(MG281)
Expression Region	1-556aa
Notes	Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity.
Tag Info	C-terminal 10xHis-tagged If you have specified tag type, please tell us and we will check if it's possible to develop.
Length/Coupling Ratio	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.