



Recombinant Human Ribitol-5-phosphate transferase FKTN (FKTN)-VLPs

Product Code	CSB-MP008709HU
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Uniprot No.	O75072
Form	Lyophilized powder
Storage Buffer	Lyophilized from PBS, 6% Trehalose, pH 7.4
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	MSRINKNVVLALLTLTSSAFLLFQLYYYYKHYLSTKNGAGLSKSKGSRIGFDSTQ WRAVKKFIMLTSNQNVPVFLIDPLILELINKNFEQVKNTSHGSTSQCKFFCVPR DFTAFALQYHLWKNEEGWFRIAENMGFQCLKIESKDPRLDGIDSLSGTEIPLHY ICKLATHAIHLVVFHERSGNYLWHGHLRLKEHIDRKFVFPFRKLQFGRYPGAFDR PELQQVTVDGLEVLIPKDPMHFVEEVPHSRFIECRYKEARAFFQQYLDDNTVE AVAFRKSARELLQLAAKTLNKLGVFPWLSSGTCLGWYRQCNIIPYSKDVDLGIF IQDYKSDIILAFQDAGLPLKHKFGKVEDSLELSFQGGKDDVKLDVFFFYEETDHM WNGGTQAKTGKKFKYLFPKFTLCWTEFVDMKVHVPCELEYIEANYGKTWKIP VKTWDWKRSPPNVQPNGIWPISSEWDEVIQLY
Research Area	Biochemicals
Source	Mammalian cell
Target Names	FKTN
Protein Names	Recombinant Human Fukutin(FKTN)
Expression Region	1-461aa
Notes	Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity.
Tag Info	C-terminal 10xHis-tagged If you have specified tag type, please tell us and we will check if it's possible to develop.
Length/Coupling Ratio	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.