



Di-Methyl-Histone H3 (Lys4) Monoclonal Antibody

Product Code	CSB-MA010418PA4me2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431/Q71DI3/P84243
Immunogen	TK(Methyl)QTARKSTGGKAPC synthetic peptide conjugate to OVA
Raised In	Mouse
Species Reactivity	Human
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	HistoneH3
Clone No.	6H12A11

Image

Western Blot

Positive WB detected in: Jurkat whole cell lysate

HistonC1 antibody at 1:1000

Secondary

Goat polyclonal to mouse IgG at 1/50000 dilution

Predicted band size: 15-25 KDa

Observed band size: 15-25 KDa

Exposure time?5min