




Recombinant Human Intestinal-type alkaline phosphatase (ALPI) (Active)

Product Code	CSB-MP001627HU
Abbreviation	Recombinant Human ALPI protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09923
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 nmol p-nitro-phenylphosphate (pNPP), in 1 minute at 37°C, pH10.0. The specific activity is > 8836.463 U/mg.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	VIPAEENPAFWNRQAAEALDAAKQLPIQKVAKNLILFLGDGLGVPTVTATRIL KGQKNGKLGPEPLAMDRFPYLALSPTYNVDRQVPDSAATATAYLCGVKANF QTIGLSAAARFNQCNTTRGNEVISVMNRAKQAGKSVGVVTTTRVQHASPAGT YAHTVNRNWYSDADMPASARQEGCQDIATQLISNMDIDVILGGGRKYMFPMG TPDPEYPADASQNGIRLDGKNLVQEWLAKHQGAWYVWNRTELMQASLDQSV THLMGLFEPGDTKYEIHRDPTLDPPLMEMTEAALRLLSRNPRGFYLFVEGGRI DHGHHEGVAYQALTEAVMFDDAIERAGQLTSEEDTLTLVTADHSHVFSFGGYT LRGSSIFGLAPSKAQDSKAYTSILYGNPGYVFNQSGVRPDVNESESGSPDYQ QQAAVPLSSETHGGEDVAVFARGPQAHLVHGVQEQSFVAHVMAFAACLEPY TACDLAPPACTTD
Source	Mammalian cell
Target Names	ALPI
Expression Region	20-503aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	55.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.