





PRLR Recombinant Monoclonal Antibody

Product Code	CSB-RA018727A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16471
Immunogen	Recombinant Human PRLR protein
Species Reactivity	Human, Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PRLR
Clone No.	1D3

Image



The Binding Activity of Human PRLR with Anti-PRLR recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human PRLR (CSB-MP018727HU1d7) at 2 µg/mL can bind Anti-PRLR recombinant antibody, the EC₅₀ is 126.8-171.9 ng/mL.



The Binding Activity of Mouse PRLR with Anti-PRLR recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Mouse PRLR at 5 µg/mL (CSB-MP018727MO) can bind Anti-PRLR recombinant antibody, the EC₅₀ is 4.021-8.706 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.