





PRLR Monoclonal Antibody

Product Code	CSB-RA018727A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16471
Immunogen	Recombinant Human PRLR protein
Species Reactivity	Human, Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Gene Names	PRLR

Image



The Binding Activity of Human PRLR with Anti-PRLR recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human PRLR (CSB-MP018727HU1d7) at 2 µg/mL can bind Anti-PRLR recombinant antibody, the EC₅₀ is 126.8-171.9 ng/mL.



The Binding Activity of Mouse PRLR with Anti-PRLR recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Mouse PRLR at 5 µg/mL (CSB-MP018727MO) can bind Anti-PRLR recombinant antibody, the EC₅₀ is 4.021-8.706 ng/mL.