




dcyD Antibody, Biotin conjugated

Product Code	CSB-PA18017D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P76316
Immunogen	Recombinant Escherichia coli D-cysteine desulfhydrase protein (2-328AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	dcyD

Image



Western Blot
Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg
All lanes: dcyD antibody at 0.25µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 36 kDa
Observed band size: 36 kDa`

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.