




cry1Bb Antibody

Product Code	CSB-PA671565XA11BRN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q45739
Immunogen	Recombinant Bacillus thuringiensis Pesticidal crystal protein Cry1Bb (390-655aa)
Species Reactivity	Bacillus thuringiensis
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus thuringiensis
Target Names	cry1Bb

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: cry1Bb antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56.7 kDa

Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.