




UL49 Antibody

Product Code	CSB-PA325857XA11HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16786
Immunogen	Recombinant Human cytomegalovirus (strain AD169) (HHV-5) (HCMV) Uncharacterized protein UL49 (354-570aa)
Species Reactivity	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Target Names	UL49

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: UL49 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50.6 kDa

Observed band size: 54 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.