



SUB1635 Antibody

Product Code	CSB-PA3334XA01FRN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B9DVS5
Immunogen	Recombinant Streptococcus uberis (strain ATCC BAA-854 / 0140J) Putative membrane anchored protein (485-495aa)
Species Reactivity	Streptococcus uberis (strain ATCC BAA-854 / 0140J)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus uberis (strain ATCC BAA-854 / 0140J)
Target Names	SUB1635
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.