



SUB1635 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA3334XA11FRN |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | B9DVS5 |
| Immunogen | Recombinant Streptococcus uberis (strain ATCC BAA-854 / 0140J) Putative membrane anchored protein (811-821aa) |
| Species Reactivity | Streptococcus uberis (strain ATCC BAA-854 / 0140J) |
| Tested Applications | ELISA, WB |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Protein G |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Streptococcus uberis (strain ATCC BAA-854 / 0140J) |
| Target Names | SUB1635 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |