




AR Antibody

Product Code	CSB-PA001975XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Recombinant Homo sapiens (Human) Androgen receptor (555-635aa)
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AR

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: AR antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35.7 kDa

Observed band size: 34 kDa