




kaiC Antibody

Product Code	CSB-PA855876XA01NHR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8YT40
Immunogen	Recombinant Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576) Circadian clock protein kinase KaiC (1-270aa)
Species Reactivity	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Target Names	kaiC

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: kaiC antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56.8 kDa

Observed band size: 69 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.