




rpaA Antibody

Product Code	CSB-PA3448XA01FUZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8Z0G7
Immunogen	Recombinant Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576) Two-component response regulator (1-247aa)
Species Reactivity	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Target Names	rpaA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: rpaA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43.3 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.