



Polyprotein Antibody

Product Code	CSB-PA3522XA11GOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C7C6C1
Immunogen	Recombinant Cassava brown streak virus Polyprotein (1430-1447aa)
Species Reactivity	Cassava brown streak virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cassava brown streak virus
Target Names	Polyprotein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.