



Polyprotein Antibody

Product Code	CSB-PA3523XA21GOB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C7C699
Immunogen	Recombinant Ugandan cassava brown streak virus Polyprotein (93-107aa)
Species Reactivity	Ugandan cassava brown streak virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ugandan cassava brown streak virus
Target Names	Polyprotein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.