



LBA Antibody

Product Code	CSB-PA355812XA21GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02238
Immunogen	Recombinant Glycine max (Soybean) (Glycine hispida) Leghemoglobin A (99-110aa)
Species Reactivity	Glycine max (Soybean) (Glycine hispida)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Target Names	LBA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.