



# LBA Antibody

<b>Product Code</b>	CSB-PA355812XA11GGV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02238
<b>Immunogen</b>	Recombinant Glycine max (Soybean) (Glycine hispida) Leghemoglobin A (10-31aa)
<b>Species Reactivity</b>	Glycine max (Soybean) (Glycine hispida)
<b>Tested Applications</b>	ELISA, WB
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Glycine max (Soybean) (Glycine hispida)
<b>Target Names</b>	LBA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.