




prfA Antibody

Product Code	CSB-PA325257XA01LPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22262
Immunogen	Recombinant <i>Listeria monocytogenes</i> serovar 1/2a (strain ATCC BAA-679 / EGD-e) Listeriolysin regulatory protein (1-237aa)
Species Reactivity	<i>Listeria monocytogenes</i> serovar 1/2a (strain ATCC BAA-679 / EGD-e)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	<i>Listeria monocytogenes</i> serovar 1/2a (strain ATCC BAA-679 / EGD-e)
Target Names	prfA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: prfA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31.4 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.