





Recombinant Human Cannabinoid receptor 1 (CNR1)-VLPs (Active)

Product Code	CSB-MP005678HU
Abbreviation	Recombinant Human CNR1 protein-VLPs (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21554
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Human CNR1 at 10 µg/mL can bind Anti-CNR1 recombinant antibody (CSB-RA005678MA01HU), the EC ₅₀ is 41.72-63.54 ng/mL.
Sequence	MKSILDGLADTTFFRTITTDLLYVGSNDIQYEDIKGDMSKLGYPQKFPLTSFRG SPFQEKMTAGDNPQLVPADQVNITEFYNKSLSSFKEENEIQCGENFMDIECF MVLNPSQQLAIAVLSLTLGTFTVLENLLVLCVILHSRSLRCRPSYHFIGSLAVAD LLGSVIFVYSFIDFHV FHRKDSRNVFLFKLGGVTASFTASVGSFLTLTAIDRYISIH RPLAYKRIVTRPKAVVAFCLMWTIAIVIAVLP LLGWNCEKLQSVCSDFPHIDETY LMFWIGVTSVLLLFIVYAYMYILWKAHSHAVRMIQRGTQKSIIIHTSEDGKVQVT RPDQARMDIRLAKTLVLILVVLICWGPELLAIMVYDVF GKMNKLIKTVFAFCMLC LLNSTVNPIIYALRSKDLRHAFRSMFPSCEGTAQPLDNSMGDSDCLHKHANNA ASVHRAAESCIKSTVKIAKVTMSVSTDTSAEAL
Source	Mammalian cell
Target Names	CNR1
Expression Region	1-472aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions)
Mol. Weight	54.6 kDa
Protein Length	Full Length
Image	



CSB-MP005678HU is detected by Mouse anti-6*His monoclonal antibody.



Activity
 Measured by its binding ability in a functional ELISA. Immobilized Human CNR1 at 10 µg/ml can bind Anti-CNR1 recombinant antibody (CSB-RA005678MA01HU), the EC₅₀ is 41.72-63.54 ng/mL.

Endotoxin	Less than 1.0 EU/ug as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.