



Recombinant Severe acute respiratory syndrome coronavirus 2 Spike glycoprotein (S) (G75V,T76I,SYLTPGD247-253del,L452Q,F490S,D614G,T859N), partial

Product Code	CSB-MP3324GMY(M7)
Abbreviation	Recombinant SARS-CoV-2 S protein (Mutant type), partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0DTC2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VNLTTRTQLPPAYTNSFTRGVYYPDKVFRSSVLHSTQDLFLPFFSNVTWFWHAIH VSGTNVIKRFDNPVLPFNDGVYFASTEKSNIIRGWIFGTTLDSTQSLIVNNAT NVVIKVCEFQFCNDPFLGVYYHKNNKSWMESEFRVYSSANNCTFEYVSQPFL MDLEGKQGNFKNLREFVFKNIDGYFKIYSKHTPINLVRDLPQGFSALEPLVDLPI GINITRFQTLALHRSSSGWTAGAAAYVGYLQPRTFLLKYNENGTITDAVDCA LDPLSETKCTLKSFTVEKGIYQTSNFRVQPTESIVRFPNITNLCPFGEVFNATRF ASVYAWNRRKRISNCVADYSVLNSASFSTFKCYGVSPTKLNDLCFTNVYADSF VIRGDEVQRQIAPGQTGKIADYNYKLPDDFTGCVIAWNSNNLDSKVGGNYNYQY RLFRKSNLKPFERDISTEIQAGSTPCNGVEGFNCYSPLQSYGFQPTNGVGYQ PYRVVVLSEFELLHAPATVCGPKKSTNLVKNKCVNFNENGLTGTGVLTESNKKF LPFQQFGRDIADTTDAVRDPQTLEILDITPCSFGGVSVITPGTNTSNQVAVLYQ GVNCTEVPVAIHADQLTPTWRVYSTGSNVFQTRAGCLIGAEHVNNSYECDIPI GAGICASYQTQTNSPRRAR
Research Area	Cancer
Source	Mammalian cell
Target Names	S
Expression Region	16-685aa(G75V,T76I,SYLTPGD247-253del,L452Q,F490S,D614G,T859N)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal 10xHis-tagged and C-terminal Flag-tagged

Mol. Weight 78.9 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.