



Recombinant Human Homeobox protein Hox-A4 (HOXA4)

Product Code	CSB-EP010654HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q00056
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTMSSFLINS NYIEPKFPPF EEY AQHSGSG GADGGPGGGP GYQQPPAPPT QHLPLQQPQL PHAGGGREPT ASYYAPRTAR EPAYPAAALY PAHGAADTAY PYGYRGGASP GRPPQPEQPP AQAKGPAHGL HASHVLQPQL PPPLQPRAVP PAAPRRCEAA PATPGVPAGG SAPACPLLLA DKSPGLGK GK EPVVYPWMKK IHVSAVNPSY NGGEPKRSRT AYTRQQVLEL EKEFHFNRYL TRRRRIEIAH TLCLSERQVK IWFQNRMKW KKDHLKLPNTK MRSSNSASAS AGPPGKAQTQ SPHLHPHPHP STSTPVPSSI
Source	E.coli
Target Names	HOXA4
Protein Names	Recommended name: Homeobox protein Hox-A4 Alternative name(s): Homeobox protein Hox-1.4 Homeobox protein Hox-1D
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.