



DEFA6 Antibody

Product Code	CSB-PA980801GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01524
Immunogen	Fusion protein of Human DEFA6
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC;ELISA: 1:2000-1:5000, WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DEFA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.