




Recombinant Human RING finger protein 207 (RNF207)

Product Code	CSB-EP019870HU
Abbreviation	Recombinant Human RNF207 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZRF8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSGAIFGPLEGPSSLDAPSIHPLVCPLCHVQYERPCLLDCFHDFCAGCLRGRA TDGRLTCPLCQHQTIVLKGPSGLPPVDRLLQFLVDSSGDGVEAVRCANCDLEC SEQDVETTYFCNTCGQPLCARCRDETHRARMFARHDIVALGQRSRDVPQKCT LHAEPYLLFSTDKLLLCIRCFRDMQKESRAHCVDLESAYVQGGERLEQAVLA VKALQTATREAIALLQAMVEEVVRHSAEEEDAIHALFGSMQDRLAERKALLLQA VQSQYEEKDKAFKEQLSHLATLLPTLQVHLVICSSFLSLANKAEFLDLGYELME RLQGIVTRPHHLRPIQSSKIASDHRAEFARCLEPLLLLGP RRVAAAASGANTLA GGLGPKALTGPHCPSPVGKMSGSPVQKPTLHRISISTKVLLAEGENTPFAEHC RHYEDSYRHLQAEMQSLKDQVQELHRDLTKHHSKAEIMGDVLHKSLQLDVQ IASEHASLEGMRVVFQEIWEEAYQRVANEQEIYEAQLHDLLQLRQENAYLTTIT KQITPYVRSIAKVKERLEPRFQAPVDEQSESLQNTHTDDSRNNAASARNNPGSV PEKREKTSEPKGNSWAPNGLSEEPLLKNMDHHRKQKNGGDVPTWREHPT
Research Area	Others
Source	E.coli
Target Names	RNF207
Expression Region	1-634aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	78.3 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.