



Recombinant Human Signal transducer and activator of transcription 3 (STAT3), Biotinylated


Product Code	CSB-EP022812HU(A4)-B
Abbreviation	Recombinant Human STAT3 protein, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40763
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	<p>MAQWNQLQQLDTRYLEQLHQLYSDSFPMELRQFLAPWIESQDWAYAASKES HATLVFHNLLGEIDQQYSRFLQESNVLYQHNLRRIKQFLQSRYLEKPMIARIV ARCLWEESRLLQTAATAAQGGQANHPTAAVVTEKQQMLEQHLQDVRKRVQ DLEQKMKVVENLQDDDFDFNYKTLKSQGDMQDLNGNNQSVTRQKMQQLEQM LTALDQMRRSIVSELAGLLSAMEYVQKTLTDEELADWKRRQQIACIGGPPNICL DRLENWITSLAESQLQTRQQIKKLEELQQKVSYKGDPIVQHRPMLLEERIVELFR NLMKSAFVVERQPCMPMHPDRPLVIKTGVQFTTKVRLLVKFPELNYQLKIKVCI DKDSGDVAALRGSRKFNILGTNTKVMNMEESNNGSLSAEFKHLTLREQRCGN GGRANCDASLIVTEELHLITFETEYVHQGLKIDLETHSLPVVVISNICQMPNAWA SILWYNMLTNNPKNVNFFTKPPIGTWDQVAEVLWSWQFSSTTKRGLSIEQLTTLA EKLLGPGVNYSGCQITWAKFCKENMAGKGFVWVWLDNIIDLVKKYILALWNE GYIMGFISKERERAILSTKPPGTFLLRFSESSKEGGVTFTWVEKDISGKTQIQSV EPYTKQQLNNMSFAEIIIMGYKIMDATNILVSPLVYLYPDIPKEEAFGKYCRPESQ EHPEADPGSAAPYLKTKFICVTPPTCSNTIDLPMSPRTLDSLMLQFGNNGEGAE PSAGGQFESLTFDMELTSECATSPM</p>
Research Area	Immunology
Source	E.coli
Target Names	STAT3
Expression Region	1-770aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-Avi-tagged



Mol. Weight 135.8 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.