




Recombinant Human Band 4.1-like protein 5 (EPB41L5)

Product Code	CSB-EP872549HUb0
Abbreviation	Recombinant Human EPB41L5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCM4
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Band 4.1-like protein 5(EPB41L5)
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MLSFFRRTLGRSSMRKHAERLREAQRAATHIPAAGDSKSIITCRVSLLDGTD VSVDLPPKAKGQELFDQIMYHLDLIESDYFGLRFMDSAQVAHWLDGTSIKKQ VKIGSPYCLHLRVKFSSEPNNLREELTRYLFVLQKQDILSGKLDPCPFDTA VQ LAAYNLQAE LGDYDLAEHSPELVSEFRFVPIQTEEMELAI FEKWKEYRGQTPA QAETNYLNKAKWLEMYGVDMHVVKARDGNDYSLGLTPTGVLVFEGDTKIGLF FWPKITRLDFKKNKLT LVVVEDDDQGKEQEHTFVFRLDHPKACKHLWKCAVE HHAFFRLRGPVQKSSHRSGFIRLGSFRYSGKTEYQT TKTNKARRSTSFERR PSKRYSRRTLQMKACATKPEELSVHNNVSTQSNGSQAWGMRSALPVSPSIS SAPVPVEIENLPQSPGTDQHDKCIPLNIDLLNSPDLLEATIGDVIGASDTMETS QALNDVNVATRLPGLGEPEVEYETLKDTSEK LKQLEMENSPLLSPRSNIDVNI N SQEEVVKLTEKCLNNVIESPGLNVMRVPDFKSNILKAQVEAVHKVTKEDSLLS HKNANVQDAATNSAVLNENNVP LPKESLETMLITPADSGSVLKEATDELDALL ASLTENLIDHTVAPQVSSTSMITPRWIVPQSGAMSNGLAGCEMLLTGKEGHGN KDGISLISPPAPFLVDAVTSSGPILAE EAVLKQKCLLTTEL
Research Area	Others
Source	E.coli
Target Names	EPB41L5
Expression Region	1-733aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	85.5 kDa


Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.