




A29L Antibody

Product Code	CSB-PA46839LA01MPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q77HM6
Immunogen	Recombinant Monkeypox virus (strain Zaire-96-I-16) (MPX) A29L protein
Raised In	Rabbit
Species Reactivity	Monkeypox virus (strain Zaire-96-I-16)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Monkeypox virus (strain Zaire-96-I-16) (MPX)

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: A29L antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 13 kDa

Observed band size: 13 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.