




Recombinant Pyrococcus horikoshii Glycerate 2-kinase (gck)

Product Code	CSB-EP530036FHX
Abbreviation	Recombinant Pyrococcus horikoshii gck protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O58231
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Pyrococcus horikoshii Glycerate 2-kinase(gck)
Immunogen Species	Pyrococcus horikoshii (strain ATCC 700860 / DSM 12428 / JCM 9974 / NBRC 100139 / OT-3)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MIAMDIREIGLRLVGEAIKAADPYRAVLNAVKV/SDDKIIVQGKEFEIKGKVYVIAL GKAACEMARAIEDILDVEDGVAVTKYGYGKELKRIKVEAGHIPDEKSILGAKE ALSILNRARENDIVFILISGGGSALFELPEEGISLEDLKLTTDLLLKSGAKIHEINTV RKHISKVKGGKLA KMIKGTGIVLIISDVVGDNLEAIASGPTVKDPTTFEDAKRILE LYDIWEKVPESVRLHIERGLRGEVEETLKEDLPNVHNFLIASNSISCEAIAREAQ RLGFKAYIMTTTLEGEAKDAGLFIGSIVQEIAERGRPFEPVVLVFGGETTVTIE GKGGKGGPNQEIALSATRKISDLEALIVAFD TDGTDGPTDAAGGIVDGT TYKKL REKGIDVEKVLKEHNSYEALKKVGGLLFTGPTGTNVNSIVIAIVTSKRGRT
Research Area	Others
Source	E.coli
Target Names	gck
Expression Region	1-440aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	54.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.