




Recombinant Escherichia coli O157:H7 Dihydrolipoyllysine-residue succinyltransferase component of 2-oxoglutarate dehydrogenase complex (sucB)

Product Code	CSB-EP006950EODb0
Abbreviation	Recombinant E.coli O157:H7 sucB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AFG7
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Escherichia coli O157:H7 Dihydrolipoyllysine-residue succinyltransferase component of 2-oxoglutarate dehydrogenase complex(sucB)
Immunogen Species	Escherichia coli O157:H7
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SSVDILVDPDPESVADATVATWHKKPGDAVVRDEVLVEIETDKVVLEVPASAD GILDAVLEDEGTTVTSRQILGRLREGNSAGKETSASKEEKASTPAQRQQASLE EQNNDALSPAIRLLAEHNLDA SAIKGTGVGGRLTREDVEKHLAKAPAKESAP AAAAPAAQPALAARSEKRVPMTLRKRVAERLLEAKNSTAMLTTFNEVNMKPI MDLRKQYGEAFEKRHGIRLGFMSFYVKAVVEALKRYPEVNASIDGDDVVYHN YFDVSMVSTPRGLVTPVLRDVDTLGMADIEKKIKELAVKGRDGKLTVEDLTG GNFTITNGGVFGSLMSTPIINPPQSAILGMHAIKDRPMAVNGQVEILPMMYLAL SYDHRLIDGRESVGFLVTIKELLEDPTRLLLDV
Research Area	Cardiovascular
Source	E.coli
Target Names	sucB
Expression Region	2-405aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	47.5 kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.