




# Recombinant Escherichia coli O157:H7 Dihydrolipoyllysine-residue succinyltransferase component of 2-oxoglutarate dehydrogenase complex (sucB)

<b>Product Code</b>	CSB-EP006950EOD
<b>Abbreviation</b>	Recombinant E.coli O157:H7 sucB protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AFG7
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Escherichia coli O157:H7 Dihydrolipoyllysine-residue succinyltransferase component of 2-oxoglutarate dehydrogenase complex(sucB)
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	SSVDILV PDLPE SVADATVATWHKKPGDAVVRDEV LVEIETDKVVLEVPASAD GILDAVLEDEGTTVTSRQILGRLREGNSAGKETS AKSEEKASTPAQRQQASLE EQNNDALSPAIRLLAEHNLDASA IKGTVGGRLTREDVEKHLAKAPAKESAP AAAAPAAQPALAA RSEKRVPMTRLRKRVAERLLEAKNSTAMLTTFNEVNMKPI MDLRKQYGEAF EKRHGIRLGFMSFYVKAVVEALKRYPEVNASIDGDDVVYHN YFDVSM AVSTPRGLVTPVLRD VDTLGMADIEKKIKELAVKGRDGKLTVEDLTG GNFTITNGGVFGSLMSTPIINPPQSAILGMHAIKDRPMAVNGQVEILPMMYLAL SYDHRLIDGRESVGFLVTIKELLEDPTRLLLLDV
<b>Research Area</b>	Cardiovascular
<b>Source</b>	E.coli
<b>Target Names</b>	sucB
<b>Expression Region</b>	2-405aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	49.9 kDa

**Protein Length**

Full Length of Mature Protein

**Image**

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.