




Recombinant Escherichia coli Poly-beta-1,6-N-acetyl-D-glucosamine N-deacetylase (pgaB)

Product Code	CSB-EP304310ENV
Abbreviation	Recombinant E.coli pgaB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P75906
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Escherichia coli Poly-beta-1,6-N-acetyl-D-glucosamine N-deacetylase(pgaB)
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	CISQSRTSFIPPQDRESLLAEQPWPHNGFVAISWHNVEDEAADQRFMSVRTS ALREQFAWLRENGYQPVSIAQIREAHRGKPLPEKAVVLTFFDDGYQSFYTRVF PILQAFQWPVAVWAPVGSWVDTPADKQVKFGDELVDREYFATWQQVREVAR RLVELASHTWNSHYGIQANATGSLLPVYVNRAYFTDHARYETAAEYRERIRLD AVKMTEYLRTKVEVNPVHVFWPYGEANGIAIEELKKLG YDMFFTLSEGLANAS QLDSIPRVLIANP SLKEFAQQIITVQEKSPQRIMHIDL D YVYDENLQQMDRNI VLIQRVKDMQISTVYLLQAFADPDGDGLVKEVWFPNRL LPMKADIFSRVAWQLR TRSGVNIYAWMPVLSWDL DPTLTRVKYLPTGEKKAQIHPEQYHRLSPFDDRVR AQVGMLYEDLAGHA AFDGILFHDDALLSDYEDASAPAITAYQQAGFSGSLSEIR QNPEQFKQWARFKSRALDFTLELSARVKAIRGPHIKTARNIFALPVIQPESEA WFAQNYADFLKSYDWTAIMAMPYLEGVAEKSADQWLIQLTNQIKNIPQAKDKS ILELQAQNWQKNGQHQAIS SQQLAHWMSLLQLNGVKNYGYYPDNFLHNQPEI DLIRPEFSTAWYPKND
Research Area	Others
Source	E.coli
Target Names	pgaB
Expression Region	21-672aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	82.6 kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.