




# Recombinant Escherichia coli O6:H1 UPF0304 protein yfbU (yfbU), Biotinylated

<b>Product Code</b>	CSB-EP364237EGX-B
<b>Abbreviation</b>	Recombinant E.coli O6:H1 yfbU protein, Biotinylated
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A8W9
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Escherichia coli O6:H1 UPF0304 protein yfbU(yfbU), Biotinylated
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MEMTNAQRLILSNQYKMMTMLDPANAERYRRLQTIIERGYGLQMRELDREFG ELKEETCRTIIDIMEMYHALHVSWSNLQDQQSIDERRVTFLGFDAATEARYLGY VRFMVNVEGRYTHFDAGTHGFNAQTPMWEKYQRMLNVWHACPRQYHLSAN EINQIINA
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	yfbU
<b>Expression Region</b>	1-164aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal MBP-tagged and C-terminal 6xHis-Avi-tagged
<b>Mol. Weight</b>	67.3 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.