




Recombinant Human Nuclear pore complex protein Nup85 (NUP85)

Product Code	CSB-EP874821HU
Abbreviation	Recombinant Human NUP85 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BW27
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Nuclear pore complex protein Nup85(NUP85)
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEELDGEPTVTLPVNSKKNQMYFDWGPGEMLVCETSFNKKKEKSEMVPSCP FIYIIRKDVVDVYSQILRKLFNESHGIFLGLQRIDEELTGKSRKSQLVRVSKNYRSV IRACMEEMHQVAIAAKDPANGRQFSSQVSILSAMELIWNLCEILFIEVAPAGPLL LHLLDWVRLHVCEVDSLADVLGSENPSKHDSFWNLVTILVLQGRLEARQML SKEADASPASAGICRIMGDLMTMPILSPGNTQTLTELELKWQHWHEECERYL QDSTFATSPHLESLLKIMLGDEAALLEQKELLSNWWYHFLVTRLLYSNPTVKPIDL HYAQAQSSLDLFLGGESSPEPLDNILLAAFEFDIHQVIKECSIALSNWWFVAHLTD LLDHCKLLQSHNLYFGSNMREFLLLEYASGLFAHPSLWQLGVDYFDYCPSELGR VSLELHIERIPLNTEQKALKVLRICEQRQMTEQVRSICKILAMKAVRNNRLGSAL SWSIRAKDAAFATLVSDRFLRDYCEGCFSDLDLIDNLGPAMMLSDRLTFLGK YREFHRMYGEKRFADAASLLLSLMTSRIAPRSFWMTLLTDALPLLEQKQVIFSA EQTYELMRCLDLTSRRPVHGSDTEQLQDDDIETTKVEMLRLSLARNLARAI REGSLEGS
Research Area	Tags & Cell Markers
Source	E.coli
Target Names	NUP85
Expression Region	1-656aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	82.5 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.