




Recombinant Escherichia coli Bifunctional (p)ppGpp synthase/hydrolase SpoT (spoT)

Product Code	CSB-EP360359ENV
Abbreviation	Recombinant E.coli spoT protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AG24
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Escherichia coli Bifunctional (p)ppGpp synthase/hydrolase SpoT(spoT)
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MYLFESLNQLIQTYLPEDQIKRLRQAYLVARDAHEGQTRSSGEPYITHPVAVAC ILAEMKLDYETLMAALLHDVIEDTPATYQDMEQLFGKSV AELVEGVSKLDKLF RDKKEAQAENFRKMIMAMVQDIRVILIKLADRTHNMRTLGS LRPDKRRRIARET LEIYSPLAHLGIHHIKTELEELGFEALYPNRYRVIKEVVKAARGNRKEMIQKILS EIEGRLQEAGIPCRVSGREKHLYSYCKMVLKEQRFHSIMDIYAFRVIVNDS DTC YRVLGQMHS LYKPRPGRVKDYAIPKANGYQSLHTSMIGPHGVPVEVQIRTED MDQMAEMGVAAHWAYKEHGETSTTAQIRAQRWMQSLLELQQSAGSSFEFIE SVKSDLFPDEIYVFTPEGRIVELPAGATPVDFAYAVHTDIGHACVGARVDRQPY PLSQPLTSGQTVEIITAPGARPNAAWLNFVSSKARAKIRQLLNKLRD DSVSL GRLLNHALGGSRKLN EIPQENIQRELD RMKLATLDDLLAEIGLGNAMSVVVAK NLQHGDASIPPATQSHGHLPIKGADGVLITFAKCCRPIPGDPIIAHVSPGKGLVI HHESCRNIRGYQKEPEKFMAVEWDKETAQEFITEIKVEMFNHQGALANL TAAI NTTTSNIQSLNTEEKDGRVYSAFIRLTARDRVHLANIMRKIRVMPDV IKVTRNRN
Research Area	Cancer
Source	E.coli
Target Names	spoT
Expression Region	1-702aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	86.8 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.