




Recombinant *Acinetobacter baumannii* Outer membrane protein Omp38 (omp38)

Product Code	CSB-EP754340ABS
Abbreviation	Recombinant <i>Acinetobacter baumannii</i> omp38 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6RYW5
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant <i>Acinetobacter baumannii</i> Outer membrane protein Omp38(omp38)
Immunogen Species	<i>Acinetobacter baumannii</i> (strain ATCC 19606 / DSM 30007 / JCM 6841 / CCUG 19606 / CIP 70.34 / NBRC 109757 / NCIMB 12457 / NCTC 12156 / 81)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ANAGVTVTPLLLGYTFQDSQHNNGGKDGNTNGPELQDDLFVGAALGIELTP WLGFEAEYNQVKGDVDGASAGAEYKQKQINGNFYVTSDLITKNYDSKIKPYVL LGAGHYKYDFDGVNRGTRGTSEEGTLGNAGVGAFWRLNDALSLRTEARATY NADEEFWNYTALAGLNVVVGGLKPAAPVVEVAPVEPTPVAPQPQELTEDLN MELRVFFDTNKSNIKDQYKPEIAKVAEKLSEYPNATARIEGHTDNTGPRKLNER LSLARANSVKSALVNEYNVDSRLSTQGFQAWDQPIADNKTKEGRAMNRRVFA TITGSRTVVVQPGQEAAAPAAAQ
Research Area	Others
Source	E.coli
Target Names	omp38
Expression Region	20-356aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	44.0 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.