




Recombinant Escherichia coli O157:H7 Outer membrane protein assembly factor BamA (bamA), partial

Product Code	CSB-EP364272EOD1
Abbreviation	Recombinant E.coli O157:H7 bamA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A942
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Escherichia coli O157:H7 Outer membrane protein assembly factor BamA(bamA),partial
Immunogen Species	Escherichia coli O157:H7
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	RNTGSFNFGIGYGTESGVSFQAGVQQDNWLGTGYAVGINGTKNDYQTYAELS VTNPYFTVDGVS LGGR LFYNDFQADDADLS DYT NKS YGTDVTLGFPINEYNSL RAGLG YVHNSLSNMQPQVAMWRYLYSMGEHPSTSDQDNSFKTDDFTFN YG WTYNKLD RGYFPTD GSRVNL TGKVTIPGSDNEY YKVTLDTATYVPIDDDHKWW VLGRTRWGYGDGLGGKEMP FYENFYAGGSSTVRGFQSN TIGPKAVYFPHQA SNYDPDYDYECATQDGAKDLCKSDDAVGGNAMAVASLEFITPTPFISDKYANS VRTSFFWDMGTVWDTNWDSSQYSGYPDYS DPSNIRMSAGIALQWMSPLGPL VFSYAQPFFKKYDGDKAEQFQFNIGKTW
Research Area	Others
Source	E.coli
Target Names	bamA
Expression Region	421-810aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	47.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.