





Recombinant Human Melanotransferrin (MFI2)

Product Code	CSB-YP013754HU
Abbreviation	Recombinant Human MELTF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08582
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Melanotransferrin(MFI2)
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE. ≥ 90% as determined by SEC-HPLC.
Sequence	GMEVRWCATSDPEQHKCGNMSEAFREAGIQPSLLCVRGTSADHCVQLIAAQE ADAITLDGGAIYEAGKEHGLKPVVGEVYDQEVGTSYYAVAVVRRSSSHVTIDTLK GVKSCHTGINRTVGWNVPVGYLVESGRLSVMGCDVLKAVSDYFGGSCVPGA GETSYSESLCRLCRGDSSGEGVCDKSPLEYYDYSGAFRCLAEGAGDVAFVK HSTVLENTDESPSRRTWTRSEEEEGECPAHEEARRTMRSSAGQAWKWAP VHRPQDESDKGEFGKRAKSRDMLG
Research Area	Others
Source	Yeast
Target Names	MELTF
Expression Region	20-302aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	58.0 kDa
Protein Length	Full Length of Mature Protein of Isoform 2
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of MELTF was greater than 90% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.