



# Recombinant

## Coxiella burnetii Hypothetical membrane spanning protein (CBU\_1863)


|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF4560DXP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Abbreviation</b>      | Recombinant Coxiella burnetii CBU_1863 protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Uniprot No.</b>       | Q83AM0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Form</b>              | Liquid or Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage Buffer</b>    | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Product Type</b>      | Recombinant<br>Coxiella burnetii Hypothetical membrane spanning protein(CBU_1863)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen Species</b> | Coxiella burnetii (strain RSA 493 / Nine Mile phase I)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purity</b>            | ≥ 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Sequence</b>          | MRNDDDTHTSTLTSSSEQVSESKIVPKDSQKPSDPNVHQNGTDVTDSAKQNAS<br>LKEVTVVTLPPDLKALEGDTHPTMPSSREVLFPNPVHKNGTPLTHSVDPNGL<br>SKDDEVTIVALESESKALEKELKKKGINYFKIGKTIANVAYTASLCLVFYKLSAPR<br>GVADGVISALVNAPSAAIFFDQFFGKVLINYNPKRVTKTHFALS LTGLAAGNIAAI<br>AGEQIAEDAVKGSSAWIFYGGLGVSLIYTFTSRTLGLPGAVGYLLNPLLGRCRP<br>FPELEQFRKDLEKNYLNPAVQRINEKTPERFIEVFYDTIKENPNKWTGYDIYKN<br>TIEKAFQVGGLIAVQHFFGLFQQLAAEGWGKISPQLENNCGLNSISAATTAFY<br>YLLVLKLPPTFRKSLEQIWILGVHNSRTLPAKIIKTALLLATLGGAGVAAYFSGSG<br>MAEETARYDQLAQNGTLPFCVNNPPLSPWEGIDGNVNELWMTIVAGVGCNFT<br>GMLFLFNPILDFFIKKSHENPWCIQSMLAVLESAPILNSNQLEKELGSWIPLIQG<br>KRPTAIQAIENKRSLVKIAENTNLPRPKFWFCGDLWKCSSNDENTRLLSETEDP<br>VPSALTA |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | CBU_1863                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Expression Region</b> | 1-603aa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tag Info</b>          | N-terminal 10xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |



**Mol. Weight** 67.7 kDa

**Protein Length** Full Length

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.