



Phospho-mFLNc (S2237) Antibody

Product Code	CSB-PA837745XA01phMO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VHX6
Immunogen	Synthesized peptide derived from Mouse FLNc around the phosphorylation site of S2237.
Species Reactivity	Mus musculus (Mouse)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Flnc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.