




TPBG Recombinant Monoclonal Antibody


Product Code	CSB-RA024093MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13641
Immunogen	Recombinant Human TPBG protein
Species Reactivity	Human, Macaca mulatta
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TPBG
Clone No.	13C3

Image



The Binding Activity of TPBG with Anti-TPBG monoclonal antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TPBG (CSB-MP024093HUB0) at 2 µg/mL can bind Anti-TPBG recombinant antibody, the EC₅₀ is 1.230-1.519 ng/mL.



The Binding Activity of Macaca fascicularis TPBG with Anti-TPBG monoclonal antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Macaca fascicularis TPBG (CSB-MP024093MOV) at 2 µg/mL can bind Anti-TPBG recombinant antibody, the EC₅₀ is 1.672-2.181 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.