




CLDN3 Recombinant Monoclonal Antibody

Product Code	CSB-RA005505MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15551
Immunogen	Recombinant Human CLDN3 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CLDN3
Clone No.	12B2

Image



The Binding Activity of CLDN3 with Anti-CLDN3 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CLDN3 (CSB-MP005505HU) at 10 µg/mL can bind Anti-CLDN3 recombinant antibody, the EC₅₀ is 23.62-34.37 ng/mL. VLPs (CSB-MP3838) is negative control.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.