




CT26/Human CCR8 Stable Cell Line

Product Code	CSB-SC004847HU3
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	P51685
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% RPMI 1640 + 10% DMSO
Resistance	Puromycin
Intended Use	RPMI 1640 + 10% FBS + 10μg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CT26 cells (green line) and transfected Human CCR8 CT26 Stable cells (red line) were stained with anti-CCR8 antibody (2μg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.