




CHOK1/Mouse CCR8 Stable Cell Line

Product Code	CSB-SC004847MO
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	P56484
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% F-12K + 10% DMSO
Resistance	Puromycin
Intended Use	F-12K + 10% FBS + 6µg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CHO-K1 cells (green line) and transfected Mouse CCR8 CHO-K1 Stable cells (red line) were stained with anti-Mouse CD198-PE antibody (2µg/1*10⁶cells) and analyzed with flow cytometry.