




CT26/Human GPC3 Stable Cell Line

Product Code	CSB-SC009705HU2
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	P51654
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% DMEM + 10% DMSO
Resistance	Puromycin
Intended Use	DMEM + 10% FBS + 10µg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CT26 cells (green line) and transfected Human GPC3 CT26 Stable cells (red line) were stained with anti-GPC3 antibody(CSB-RA009705A1HU) (2µg/1*10⁶cells), washed and then followed by APC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.