




CHOK1/Human CD69 Stable Cell Line

Product Code	CSB-SC004952HU2
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	Q07108
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% F-12K + 10% DMSO
Resistance	Puromycin
Intended Use	F-12K + 10% FBS + 10µg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CHO-K1 cells (green line) and transfected Human CD69 CHO-K1 stable cells (red line) were stained with anti-CD69 antibody (biolegend/310939/FN50) (2µg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Mouse IgG Fc antibody and analyzed with flow cytometry.