




HEK293T/Human LY6G6D Stable Cell Line

Product Code	CSB-SC013246HU2
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	O95868
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% DMEM + 10% DMSO
Resistance	Puromycin
Intended Use	DMEM + 10% FBS + 2μg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected HEK293T cells (green line) and transfected Human LY6G6D HEK293T stable cells (red line) were stained with anti-LY6G6D antibody (CSB-RA013246MA1HU) (2μg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.